

## Deha Oxygen/Hydrazine Test Kit-Comparator

### How to use:

The Comparators rely on visual assessment. The measurement result is not detected electronically. The human eye compares the measurement colour that results from adding reagents (e.g. tablets) with a comparative scale (colour disc)

Measurement-Range extendable by dilution: DEHA 0.00 - 0.50 mg/l

### Reagents:

30ml DEHA Test solution 65ml DEHA Indicator sol. (125 tests)

### Supplied with:

Black plastic carrying case with foam insert  
2 x 10ml plastic square cell  
1 x plastic stirring rod + brush

