

## Sulphite (HR) Test Kit – Tablet Count Method

### How to use:

A tablet is added to a 50ml or 100ml sample. The water sample develops a colour. More tablets are added until a colour change happens. The number of tablets used to force the colour change are entered into a simple formula to obtain the test result.

### Measurement-Range extendable by dilution:

Sulphite  
(10 - 500 mg/l Na SO)

### Reagents:

Sulphite HR No. 1 (100 tabl.)  
Sulphite HR No. 2 HR (100 tabl)

### Supplied with:

- Black plastic carrying case with foam insert
- 1 x 100ml shaker tube
- plastic stirring rod
- cleaning brush
- printed manual
- SVS2800 with 3 shaker tubes

